

## TAP-Protein prep for *in vitro* kinase assay

5/28-2015 (modified by Ikui)

### Cultures:

1. Grow O/N culture to O.D. (0.02-0.05) in 25 ml Sc-ura + 3% raffinose. Grow O/N in 30°C
2. Dilute O/N culture to O.D (0.3-0.6) in 100 ml Sc-ura + 3% Raffinose. Grow to O.D. (0.5-0.8).
  - a. Pre-induction samples: take 5 ml pellet, and freeze in -80°C
3. Add 3% galactose, grow to O.D 1.0. Take 5ml post-induction samples.
4. Spin samples in 250 ml bottles at 3000 rpm, 4°C for 10 minutes
5. Transfer pellets and wash with 10 ml cold H<sub>2</sub>O into 50 ml conical tubes
6. Spin at 3500rpm, 4°C for 6 minutes.
7. Pellet and freeze in -80°C

### Lysis and IP:

1. Pre-induction sample extraction:
  - a. Transfer samples (step 3a above) to spin cap tubes with 1ml TE and pellet.
  - b. Add 2 scoops glass beads and 300µl 1x SDS onto pellet. Fast prep twice.
  - c. Heat at 95°C for 5 min. Spin for 2 min.
2. Post-induction sample (step 4 above) extraction: Same as step 1.
3. Resuspend pellets (step 8 above) in 300ul TAP lysis buffer (+ fresh inhibitors)
  - a. Add 2 scoops of glass beads into cap tubes per sample
  - b. Fast prep (setting 6) for 60 seconds **twice**. (wait 5 minutes between the shake) Spin, then pool 6 supernatants in a tube.
  - c. Collect 20µl samples for western + add 20µl 2x sample buffer
4. Spin down lysates for 1 hr in 10,000 rpm, 4°C. Collect 20µl from supernatant for western
5. Equilibrate rabbit IgG agarose beads (Sigma) and add samples:
  - a. Prepare 1ml TAP Core + fresh inhibitors per sample
  - b. Equilibrate 40µl beads in 1ml TAP core buffer for each sample in spin cap tubes
    - Spin at 3000 rpm, discard supernatant
  - c. Add supernatant from samples into beads
  - d. Put tubes on rocker in 4°C for 2 hours
  - e. Spin down beads for 5 minutes at 1500 rpm, 4°C. Take 20µl sample from supernatant for western
  - f. Wash with 1.5 ml TAP CORE buffer, and spin for 15 min at 4°C. Repeat 5 times
6. Cleave off TAP bound to beads using TEV:
  - a. Wash with 1ml TEV cleavage buffer, spin for 10 min, RT. Remove supernatant
  - b. Add 30µl TEV buffer, 2µl TEV protease. Incubate in 30°C for 30 minutes
  - c. Spin down. Aliquot 10µl from each into eppendorf tubes. Take 1µl sample for western. Add 10µl 1x SDS Buffer
  - d. Store purified proteins in -80°C.
  - e. (if you use 3C protease from Clontech, incubate at 4 degrees ON)

7. wash beads with 1ml H<sub>2</sub>O, and add 20ul 1x sample buffer to tube with beads for western.

**Reagents:**

TAP lysis Buffer (200 ml stock, 50mM HEPES 7.4, 300mM NaCl, 0.1% NP40, 50mM Naf, 80mM BGP)


|                         |   |
|-------------------------|---|
| 20ml                    | 0.5M HEPES pH 7.4 (Sigma H3375)   |
| 12 ml                   | 5M NaCl   |
| 200µl                   | 100% NP40 (Fisher NC9375914)  |
| 20ml                    | 500mM NaF (Sigma S7920)   |
| 16ml                    | 1M BGP (Sigma G9422)  |
| 131.8ml                 | Autoclaved dH <sub>2</sub> O  |
| <i>Fresh before use</i> | Protease inhibitors, Phosphatase Inhibitors, PMSF, Aprotinin<br>(PMSF Sigma P7626, Aprotinin Sigma A6279) |

TAP Core (500 ml, 50mM HEPES 7.4, 150mM NaCl, 0.1% NP40, 50mM Naf, 80mM BGP)

|                         |  |
|-------------------------|--|
| 50ml                    | 0.5M HEPES pH 7.4 (Sigma H3375)                                |
| 15ml                    | 5M NaCl  |
| 500µl                   | 100% NP40 (Fisher NC9375914)                                   |
| 50ml                    | 500mM NaF (Sigma S7920)  |
| 40ml                    | 1M BGP (Sigma G9422)   |
| 344.5ml                 | Autoclaved dH <sub>2</sub> O                                   |
| <i>Fresh before use</i> | Aprotinin, PMSF (PMSF: Sigma P7626, Aprotinin: Sigma<br>A6279) |

TEV cleavage Buffer (200ml)

|                        |                                  |
|------------------------|----------------------------------|
| 20ml                   | 500mM HEPES pH 8.0 (Sigma H3375) |
| 6ml                    | 5M NaCl                          |
| 200µl                  | 100% NP40 (Fisher NC9375914)     |
| 200µl                  | 500mM EDTA (Fisher BP118)        |
| 20ml                   | 100% Glycerol (Fisher G33)       |
| 153.4ml                | Autoclaved dH <sub>2</sub> O     |
| 200ul (add before use) | 1M DTT (Sigma D0632)             |



DOM 638 – Mck1-TAP IP  
Probed with anti-PrA-HRP, 1:5000, 1 hr  
Exposure: 1 min, ECL regular  
Mck1-TAP expected size: 62 kDa  
07/16/2013